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RIC Rawalpindi
Microbiology Chart Report
Printed May 11, 2019 20:03 PKT

bioMérieux Customer: Micro Lab

Patient Name: _____ Patient ID: _____
Location: _____ Physician: _____
Lab ID: 0091 Isolate Number: 1

Organism Quantity: _____
Selected Organism : **Francisella tularensis**

Source: _____ Collected: _____

Comments: _____


Identification Information	Analysis Time: 9.95 hours	Status: Final
Selected Organism	92% Probability Bionumber: 0002000100001200	Francisella tularensis
ID Analysis Messages	Confirm by serological tests Highly pathogenic organism	

2	APPA	-	3	ADO	-	4	PyrA	-	5	IARL	-	7	dCEL	-	9	BGAL	-
10	H2S	-	11	BNAG	-	12	AGLTp	-	13	dGLU	-	14	GGT	+	15	OFF	-
17	BGLU	-	18	dMAL	-	19	dMAN	-	20	dMNE	-	21	BXYL	-	22	BAIap	-
23	ProA	-	26	LIP	-	27	PLE	-	29	TyrA	+	31	URE	-	32	dSOR	-
33	SAC	-	34	dTAG	-	35	dTRE	-	36	CIT	-	37	MNT	-	39	5KG	-
40	ILATk	-	41	AGLU	-	42	SUCT	-	43	NAGA	-	44	AGAL	-	45	PHOS	-
46	GlyA	+	47	ODC	-	48	LDC	-	53	IHISa	-	56	CMT	+	57	BGUR	-
58	O129R	-	59	GGAA	-	61	IMLTa	-	62	ELLM	-	64	ILATa	-			

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+= Deduced drug *= AES modified **= User modified

AES Findings	
Confidence:	Analysis not performed

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
Infectious diseases	Annual deaths (millions)
Respiratory infections	3.96
HIV/AIDS	2.77
Diarrhoeal diseases	1.80
Tuberculosis	1.56
Vaccine-preventable childhood diseases	1.12
Malaria	1.27
STDs (other than HIV)	0.18
Meningitis	0.17
Hepatitis B and C	0.16
Tropical parasitic diseases	0.13
Dengue	0.02
Other infectious diseases	1.76

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Top 10 global causes of deaths, 2016

Source: Global Health Estimates 2016: Deaths by Cause, Age, Sex, by Country and by Region, 2000-2016. Geneva, World Health Organization, 2018.

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
Routine molecular point-of-care testing for respiratory viruses in adults presenting to hospital with acute respiratory illness (ResPOC): a pragmatic, open-label, randomised controlled trial

Nathan J Brendish, Ahalya K Malachira, Lawrence Armstrong, Rebecca Houghton, Sandra Aitken, Esther Nyimbili, Sean Ewings, Patrick J Lillie, Tristan W Clark

Summary

Background Respiratory virus infection is a common cause of hospitalisation in adults. Rapid point-of-care testing (POCT) for respiratory viruses might improve clinical care by reducing unnecessary antibiotic use, shortening length of hospital stay, improving influenza detection and treatment, and rationalising isolation facility use; however, insufficient evidence exists to support its use over standard clinical care. We aimed to assess the effect of routine POCT on a broad range of clinical outcomes including antibiotic use.

Methods In this pragmatic, parallel-group, open-label, randomised controlled trial, we enrolled adults (aged ≥ 18 years) within 24 h of presenting to the emergency department or acute medical unit of a large UK hospital with acute respiratory illness or fever higher than 37.5°C (≤ 7 days duration), or both, over two winter seasons. Patients were randomly assigned (1:1), via an internet-based allocation sequence with random permuted blocks, to have a molecular POC test for respiratory viruses or routine clinical care. The primary outcome was the proportion of patients who received antibiotics while hospitalised (up to 30 days). Secondary outcomes included duration of antibiotics, proportion of patients receiving single doses or brief courses of antibiotics, length of stay, antiviral use, isolation facility use, and safety. Analysis was by modified intention to treat, excluding patients who declined intervention or were withdrawn for protocol violations. This study is registered with ISRCTN, number 90211642, and has been completed.

Findings Between Jan 15, 2015, and April 30, 2015, and between Oct 1, 2015, and April 30, 2016, we enrolled 720 patients (362 assigned to POCT and 358 to routine care). Six patients withdrew or had protocol violations. 301 (84%) of 360 patients in the POCT group received antibiotics compared with 294 (83%) of 354 controls (difference 0.6%, 95% CI -4.9 to 6.0 ; $p=0.84$). Mean duration of antibiotics did not differ between groups (7.2 days [SD 5.1] in the POCT group vs 7.7 days [4.9] in the control group; difference -0.4 , 95% CI -1.2 to 0.4 ; $p=0.32$). 50 (17%) of 301 patients treated with antibiotics in the POCT group received single doses or brief courses of antibiotics (<48 h) compared with 26 (9%) of 294 patients in the control group (difference 7.8%, 95% CI 2.5 to 13.1 ; $p=0.0047$; number needed to test=13). Mean length of stay was shorter in the POCT group (5.7 days [SD 6.3]) than in the control group (6.8 days [7.7]; difference -1.1 , 95% CI -2.2 to -0.3 ; $p=0.0443$). Appropriate antiviral treatment of influenza-positive patients was more common in the POCT group (52 [91%] of 57 patients) than in the control group (24 [65%] of 37 patients; difference 26.4%, 95% CI 9.6 to 43.2 ; $p=0.0076$; number needed to test=4). We found no differences in adverse outcomes between the groups (77 [21%] of 360 patients in the POCT group vs 88 [25%] of 354 patients in the control group; -3.5% , -9.7 to 2.7 ; $p=0.29$).

Interpretation Routine use of molecular POCT for respiratory viruses did not reduce the proportion of patients treated with antibiotics. However, the primary outcome measure failed to capture differences in antibiotic use because many patients were started on antibiotics before the results of POCT could be made available. Although POCT was not associated with a reduction in the duration of antibiotics overall, more patients in the POCT group received single doses or brief courses of antibiotics than did patients in the control group. POCT was also associated with a reduced length of stay and improved influenza detection and antiviral use, and appeared to be safe.

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Bakteriyalar

Yarı-Kantitatif

Acinetobacter calcoaceticus-baumannii complex

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
www.BioFireDx.com

Run Information

Sample ID	123-xyz-sampID	Run Date	05 Sept 2017 8:00 AM
Protocol	BAL v2.0	Serial No.	01234567
Pouch Type	Pneumo v2.0	Lot No.	123456
Controls	Passed	Operator	Jane Doe (JD)
Run Status	Completed	Instrument	FA1234

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


1-MacVane, SH and Nolte, FS. JCM, 2016Mar;54(10):2455–63.;
2- Banerjee, R, et al. Clinical Infectious Diseases. 2015;61(7):1071–80;
3- Ray et al. The Pediatric Infectious Disease Journal • Volume 35, Number 5, May 2016;
4-Pardo et al. Diagnostic Microbiology and Infectious Disease 84 (2016) 159–164

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
European Journal of Clinical Microbiology & Infectious Diseases (2020) 39:293–297
<https://doi.org/10.1007/s10096-019-03724-7>

ORIGINAL ARTICLE

Impact of a multiplex PCR assay (FilmArray®) on the management of patients with suspected central nervous system infections

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Received: 31 July 2019 / Accepted: 24 September 2019 / Published online: 12 November 2019
© Springer-Verlag GmbH Germany, part of Springer Nature 2019

Abstract

Microbiological diagnosis of central nervous system (CNS) infections is challenging due to limited access to CNS samples, overlap between meningitis and encephalitis, and the multiplicity of pathogens potentially involved. We aimed to estimate the impact of a commercial multiplex PCR assay (FilmArray® meningitis/encephalitis) on the management of patients with suspicion of meningitis or encephalitis, in terms of time to diagnosis, antimicrobial agents use, duration of hospitalization, and costs. This prospective observational study was conducted at Saint Joseph Hospital (Paris, France) from December 2016 to December 2017. All CSF samples sent to the microbiology laboratory for suspicion of meningitis and/or encephalitis, with CSF cells count > 5 cells/μL, were tested by meningitis/encephalitis multiplex PCR assay. One hundred thirty patients were included. The multiplex PCR assay was positive in 33 patients (25%). Main pathogens found were Enterovirus ($n = 12$), Varicella-Zoster virus ($n = 7$), Herpes simplex virus-2 ($n = 6$), and *Listeria monocytogenes* ($n = 3$) as main pathogens. The multiplex PCR assay reduced time to microbiological diagnosis by 3.3 ± 1.6 days and allowed an earlier discontinuation of empirical anti-infective drugs in 42 patients (32%) and an earlier hospital discharge in 23 patients (18%), with an estimated saving of 82 hospital days overall, and a management cost reduction of 26,242 € (201 €/patient). The systematic use of the FilmArray® meningitis/encephalitis multiplex PCR assay may allow earlier diagnosis, earlier discontinuation of empirical treatment, reduced duration of stay, and costs reduction.

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